

Quick Start Protocol

GelRed® & GelGreen® Nucleic Acid Gel Stains

Catalog nos. 41002, 41003, 41004, 41005

Precast agarose gels

1. Add GelRed® or GelGreen® to molten agarose at 1X final concentration.

For example, add 5 μ L of 10,000X GelRed® or GelGreen® to 50 mL agarose.

2. Cast gel and run samples.

Recommended loading is 50-200 ng DNA or ladder per lane, or 2-5 μ L PCR product. If DNA concentration is unknown, run 1/2 to 1/3 the amount you would normally load on an EtBr gel. Overloading DNA can lead to poor band migration. If you need to load more DNA, use the post-staining protocol (other side).

3. Image gel.

GelRed®: use UV light box with EtBr filter.

GelGreen®: use UV light box with SYBR® Green filter, or blue light illuminator such as a Dark Reader®.

Visit www.biotium.com for complete protocols, troubleshooting tips, and answers to frequently asked questions (FAQs).

Rev: November 9, 2018



Post-staining agarose or acrylamide gels

1. **Run samples on gel with no DNA gel stain added.**
2. **Dilute GelRed® or GelGreen® in water at 3X final concentration.**
For example, add 15 μ L 10,000X GelRed® or GelGreen® to 50 mL water.
3. **Place gel in clean container with enough 3X gel stain to cover gel and incubate with rocking for 30 min.**
Bands may begin to be detectable after 5 min.
4. **Image gel.**
GelRed®: use UV light box with EtBr filter.
GelGreen®: use UV light box with SYBR® Green filter, or blue light illuminator such as a Dark Reader®.

Visit www.biotium.com/support for complete protocols, troubleshooting tips, and answers to frequently asked questions (FAQs).

For research use only. GelRed® and GelGreen® are covered by US and international patents. Dark Reader is a registered trademark of Clare Chemical; SYBR is a registered trademark of Thermo Fisher Scientific.