## Lucigen® Simplifying Genomics

## LavaLAMP™ DNA Component Kits

(CAT. NO. 30076-1, 30077-1)

- 1. Thaw all kit components on ice and set up reactions on ice.
- 2. Mix each component thoroughly before use by vortexing for 3 10 sec. Centrifuge briefly to collect contents.
- 3. Prepare reaction mix(es) in the order listed below (Table 1).

## Notes:

- The reaction conditions recommended (Table 1) are for the use of Green Fluorescent Dye included with Cat. No. 30077-1 and 1  $\mu$ L of Target DNA Sample. Adjust the volume of nuclease-free H<sub>2</sub>O when using other dye or target DNA sample amounts.
- Table 1 provides volumes for a single reaction, if multiple reactions are required increase volumes proportionately.

  Prepare enough reaction mix cocktail(s) for the number of amplification reactions being performed plus an additional 10% to accommodate slight pipetting errors.

Table 1. Positive Control, No Target Control and Experimental Reaction Setup (suggested)

· · · · · · · · · · · · · · · · · · ·		. , 50 /	
Component	Positive Control	No Target Control (NTC)	Experimental
	Amount (μL)	Amount (μL)	Amount (μL)
Nuclease-free H <sub>2</sub> O	13.7	13.7	13.7
10X LavaLAMP™ DNA Buffer	2.5	2.5	2.5
LavaLAMP™ DNA Enzyme	1.0	1.0	1.0
dNTP mix, 25 mM	0.8	0.8	0.8
Magnesium Sulfate, 100 mM	2.5	2.5	2.5
Green Fluorescent Dye (optional)	1.0	1.0	1.0
Target-Specific Primer Mix, 10X		2.5	2.5
DNA Positive Control LAMP Primer Mix	2.5		
Total	24.0	24.0	24.0

- 4. Mix the reagents completely by pipetting.
- 5. If more than one reaction is being run, dispense 24  $\mu$ L of the reaction mix for each reaction into PCR tubes or a 96-well PCR plate.

**Note**: To minimize cross-contamination, perform steps 6-8 in an area separate from the area used to assemble the reaction mix.

- 6. Add 1  $\mu$ L of Target DNA or Positive Control DNA to the appropriate reaction tubes or wells. Add 1  $\mu$ L of nuclease-free water to the NTC reaction tubes or wells. Mix completely by pipetting.
- 7. Cap tubes or seal plate wells. Centrifuge briefly to collect contents prior.
- 8. Incubate the reactions as follows:

Step	Temperature	Time
1. Amplification	<b>Experimental and NTC</b> : 68°C – 74°C <b>Positive Control</b> : 74°C	30 - 60 minutes
2. Hold (Optional)	4°C	∞

- 9. Immediately stop the amplification reactions by inactivation, using one of the three methods below.
  - a. Place on ice or at 4°C.
  - b. Add gel loading dye to yield 10 mM EDTA (final concentration.)
  - c. Heat-kill in a thermocycler or heat block at 95°C for 5 minutes
- 10. Detect amplified product using your method of choice.