


# Optimize Your Transfections!


5 tips to get the best results in your experiments

## Get Your Ratios Right

It's important to test multiple **Reagent:DNA ratios**. Finding the optimal ratio for your cells will help unlock the highest transfection efficiency and decrease toxicity.



## The Complexity of Complexes



DNA-Reagent complex formation takes time. Allow your transfection reagent and DNA mixture to incubate for 15-30 minutes\* at room temperature before adding the mix to your cells. Transfection efficiency will decrease if complex formation time exceeds one hour.

\*Incubation time varies by reagent, reference product protocol

## DNA Do's and Don'ts

**Do:** Use endotoxin-free, **high quality DNA** for all transfections. Contaminants like proteins, carbohydrates and lipids can negatively affect transfection efficiency and gene expression levels.

**Don't:** Assume all cell types require the same dose of DNA. In some suspension cell lines, such as THP-1, transfecting higher amounts of DNA leads to higher transgene expression. Be sure to keep DNA:Reagent ratio constant when adjusting DNA dose.

## Divide and Conquer

Use **happy, healthy, and actively dividing cells** to maximize your transfection efficiency. Mirus recommends plating your cells the night before transfection at a density that will promote cell division, and obtain **75-90% confluency** at the time of transfection.



Visit [mirusbio.com/optimize](https://mirusbio.com/optimize)

• Full optimization protocol • More transfection tips  
• [techsupport@mirusbio.com](mailto:techsupport@mirusbio.com) | 844.647.8724