

INSTRUCTION MANUAL

Zymoprep™ Yeast Plasmid Miniprep I Catalog No. D2001

Highlights

- Simple procedures for plasmid rescue from yeast.
- Ideal for low copy and hard to isolate plasmids.
- For isolation of plasmid DNA for downstream applications such as PCR, transformation, hybridization, etc.

Contents

| Product Contents | 1 |
|----------------------|---|
| Specifications | 1 |
| Product Description | 2 |
| Protocol | |
| Ordering Information | 4 |
| Related Products | 5 |

Satisfaction of all Zymo Research products is guaranteed. If you should be dissatisfied with this product please call 1-888-882-9682. For Research Use Only

Product Contents

| Zymoprep™ Yeast Plasmid Miniprep I (Kit Size) | D2001 (100 preps) | Storage Temperature |
|--|--------------------------|------------------------|
| Solution 1, Digestion Buffer | 15 ml | RT |
| Solution 2, Lysis Buffer | 15 ml | RT |
| Solution 3, Neutralizing Buffer | 15 ml | RT |
| Zymolyase™ and Storage Buffer* | 1,000 units | -20°C |
| Instruction Manual | 1 | - |

Note - Integrity of kit components is guaranteed for up to one year at 4°C from date of purchase. Reagents are routinely tested on a lot-to-lot basis to ensure they provide the highest performance and reliability.

Specifications

- Sample Sources S. cerevisiae, C. albicans and S. pombe, and other fungi species sensitive to yeast lytic enzymatic digestion (Zymolyase™).
- Format Isopropanol precipitation.
- Plasmid Size DNA up to 23 kb.
- Equipment Needed Incubator shaker, microcentrifuge

Note - TM Trademarks of Zymo Research Corporation. This product is for research use only and should only be used by trained professionals. This product is not intended for use in diagnostic procedures. Some reagents included with this kit are irritants. Wear

^{*}The Zymolyase™ is stable as shipped. Add 200 µl of supplied Storage Buffer to each Zymolyase™ tube prior to use. The final concentration of Zymolyase™ after the addition of the Storage Buffer is 5 units/µl.

protective gloves and eye protection. Follow the safety guidelines and rules enacted by your research institution or facility.

Product Description

The **Zymoprep**TM is a simple and efficient yeast plasmid miniprep that is based on the *E. coli* alkaline lysis method but using **Zymolyase**TM as the first solution. There is no need for glass beads, phenol, or vortexing. Instead, plasmid DNA is reliably recovered from yeast cells whether colonies, patches on plates, or liquid cultures are sampled. Plasmid yields are typically between 0.01-0.3 ng for most 2 μ based plasmids from 1.5 ml overnight cultures. This kit also works well with low copy number yeast plasmids. Recovered plasmid DNA is in TE buffer and can be used for *E. coli* transformation, Western blotting, PCR, etc.

For **Assistance**, please contact Zymo Research Technical Support at 1-888-882-9682 or e-mail tech@zymoresearch.com.

Reagent Preparation:

✓ Add 200 µl of the supplied Storage Buffer to the lyophilized Zymolyase™. Mix to dissolve the enzyme completely and spin down briefly using a microcentrifuge. Store the reconstituted Zymolyase™ at -20°C.

Standard Protocol

Grow yeast cells at 30°C in YPD broth or selective medium. Unless stated otherwise, the following steps in the procedure are performed at room temperature.

- 1. Aliquot 0.5-1.0 ml of the full-grown yeast cells into 1.5 ml microcentrifuge tubes and spin down the cells at 600 x g for 2 minutes. Discard the supernatant.
- 2. Add 150 µl Solution 1 to each pellet.
- 3. Add 2 µl of **Zymolyase**™ to each tube. Resuspend the pellet by flicking the tube with your finger or vortexing.

Note: For multiple sample processes, add 13 µl Zymolyase for each ml of Solution 1 to make a Solution 1-enzyme mixture. Use 150 µl of this mixture to re-suspend the pellet for each sample.

- 4. Incubate at 37°C for 15-60 minutes (15 minutes is the minimal incubation time, although longer incubation is optional).
- 5. Add 150 µl **Solution 2** to each tube. Mix well.
- 6. Add 150 µl **Solution 3** to each tube. Mix well.
- 7. Centrifuge at maximum speed for 2 minutes.
- 8. Transfer the supernatant to new tubes. Add 400 µl isopropanol (2-propanol) to each tube. Mix well.
- 9. Centrifuge at maximum speed for 8 minutes. Aspirate supernatant. Spin briefly again and remove any residual supernatant.
- 10. Resuspend the plasmid pellet in 35 µl TE buffer. It is not necessary to dry the pellet before adding the TE. Sometimes the pellet requires repeated pipetting to be completely dissolved.

| Page 3 | |
|--------|--|
| | Use 3-5 µl of the plasmid DNA for <i>E. coli</i> transformation experiments. |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |

Protocol for use with colonies or patches

- Use toothpick or pipette tip to pick roughly 5-15 µl volume of yeast colonies or patches from plates and dispense into 150 µl of Solution 1-enzyme mixture (add 13 µl Zymolyase™ to each ml of Solution 1 to make Solution 1-enzyme mixture).
- 2. Incubate at 37°C for 15-60 minutes (15 minutes is the minimal incubation time, although longer incubation is optional).
- 3. Add 150 µl Solution 2 to each tube. Mix well.
- 4. Add 150 µl **Solution 3** to each tube. Mix well.
- 5. Centrifuge at maximum speed for 2 minutes.
- 6. Transfer supernatant to new tubes. Add 400 μl isopropanol (2-propanol) to each tube. Mix well.
- 7. Centrifuge at maximum speed for 8 minutes. Aspirate supernatant. Spin briefly and remove any residual supernatant.
- Resuspend the plasmid pellet in 35 μl TE buffer. It is not necessary to dry the
 pellet before adding TE. Sometimes the pellet needs to be pipette for complete
 dissolving.

Use 3-5 µl of the plasmid DNA for *E. coli* transformation experiments.

Ordering Information

| Product Description | Catalog No. | Kit Size |
|-------------------------------------|-------------|------------|
| Zymoprep™ Yeast Plasmid Miniprep I | D2001 | 100 preps. |
| Zymoprep™ Yeast Plasmid Miniprep II | D2004 | 50 preps. |

| For Individual Sale | Catalog No. | Amount |
|-------------------------------|-------------|------------|
| Zymolyase™ and Storage Buffer | E1004 | 1000 units |
| (lyophilized) | E1005 | 2000 units |

Other Popular Yeast Purification Products from Zymo Research

| Product | Format | Kit Size | Cat No. |
|--|---|--|--------------------------------------|
| | Yeast Growth & Transformation | | |
| Frozen-EZ Yeast Transformation II Kit™ | Transformation efficiency 10 ⁵ -10 ⁶ CFU/µg | 120 rxns | T2001 |
| YPD Plus™ | Increases yeast transformation efficiency >50% | 50 ml 100 ml | Y1003-50 Y1003-100 |
| | Yeast Specialty Products | · | |
| Yeast Protein Kit™ | Efficient lysis of yeast for downstream protein and DNA analyses | 200 preps. | Y1002 |
| 5-Fluorootic Acid (5-FOA) | Yeast genetic counter selection | 1 g. 5 g. 250 ml (2X SC) 10 ml (100X) | F9001-1 F9001-5 F9002 F9003 |
| α-Factor Mating Pheromone | Optimized for yeast mating induction | 240 ul | Y1001 |
| Zymolyase-Yeast Lytic Enzyme | Efficient digestion of yeast and fungal walls | 1,000 U 2,000 U | E1004 E1005 |
| | Yeast DNA/RNA Purification | | |
| ZymoPrep™ Yeast Plasmid Miniprep I | Isopropanol precipitation Format, elution ≥ 35 μI | 100 preps. | D2001 |
| ZymoPrep™ Yeast Plasmid Miniprep II | Spin Column Format (up to 5 µg/prep.) | 50 preps. | D2004 |
| YeaStar™ Genomic DNA Kit | Spin Column Format (up to 20 μg/prep.) | 40 preps. | D2002 |
| ZR Soil Microbe DNA MiniPrep™ | Bead Bashing, Spin Column Format (up to 25 µg/prep.) | 50 preps. | D6001 |
| ZR Fungal/Bacterial DNA MiniPrep™ | Bead Bashing, Spin Column Format (up to 25 µg/prep.) | 50 preps. | D6005 |
| ZR Fungal/Bacterial RNA MiniPrep™ | Bead Bashing, Spin Column Format (up to 25 µg/prep.) | 50 preps. | R2014 |
| YeaStar™ Genomic RNA Kit | Spin Column Format (up to 25 µg/prep.) | 40 preps. | R1002 |