


# Application Note

**Generating  
growth factor gradients  
with PODS<sup>®</sup> crystals**



# Achieving growth factor gradients utilizing PODS<sup>®</sup> BMP-2

Data Courtesy of Takaki Shima, HI-LEX Corporation  
Takarazuka, Hyōgo, Japan (Issued March, 10<sup>th</sup> 2019)

## Introduction to PODS<sup>®</sup>

### The challenge with soluble growth factors

Many proteins, especially growth factors and cytokines, when used as a reagent, degrade quickly, rapidly losing their bioactivity. This fragility hampers research and significantly limits the therapeutic potential of proteins.

### Protein Micro-depots


Development of a technology that can continuously replenish active protein from a local, microscopic store, has been a significant challenge, but one that could transform the fields of cell culture and medicine by allowing greater control over the growth of cells.

### Introducing PODS<sup>®</sup>

PODS<sup>®</sup> technology has made the goal of a micro-depot for proteins a reality. PODS<sup>®</sup> is a sustained release system which continuously replenishes proteins from millions of local microscopic stores which can be placed next to (or at a distance from) cells, either randomly or in precise locations. Just like cells, these micro-depots release a steady stream of bioactive protein. This protein can be limited to local surroundings or dispersed more widely, or made to form a gradient.

### How does it work?

At the heart of PODS<sup>®</sup> is an extraordinary polyhedrin protein. This specific polyhedrin protein has the unique ability to encase cargo proteins within perfect, transparent, cubic, micro-sized crystals, much smaller than the cells. These protein crystals form admixtures of the polyhedrin and cargo proteins which slowly degrade releasing the biologically active cargo protein.




### How can PODS<sup>®</sup> help my research?

PODS<sup>®</sup> are tough and will withstand physical and chemical stress, so you can handle them with ease. PODS<sup>®</sup> can be made to release intact cargo protein over days, weeks or even months. Using PODS<sup>®</sup> you can readily create a steady-state protein environment in microscopic detail wherever you want, tailored exactly to your requirements. This is the power of PODS<sup>®</sup>. PODS<sup>®</sup> proteins are now available for many growth factors and cytokines and are already being used in many leading world-class research labs. PODS<sup>®</sup> protein applications include:

- Micropatterning
- Physiological, stable gradient formation
- Bioinks for 3D printing
- Microcarriers
- Functionalizing scaffolds
- Microfluidics (lab on a chip)
- Improved and simplified stem cell culture
- Therapeutic protein delivery

## Methods




**Culture method:** PODS<sup>®</sup> BMP-2 and PODS<sup>®</sup> Empty crystals were spotted as six dense discs into wells of a 24-well plate and dried on (black squares). ATDC5 cells, a chondrogenic cell line, were then seeded across the entire well and cultured for 9 days. **NOTE:** a single application of PODS<sup>®</sup> crystals was used during the culture period.

**Staining:** On Day 9, cells were fixed in methanol. Subsequently, cells were either stained with Alizarin red for calcification, or stained with Alcian blue for extracellular matrix (ECM) production.

## Results

**A**

PODS <sup>®</sup> BMP-2	PODS <sup>®</sup> Empty	Cells only
----------------------------	----------------------------	---------------




**(A) Alizarin red staining of ATDC5 cells**

ATDC5 cells stained with Alizarin red to indicate osteogenic differentiation (calcification). PODS<sup>®</sup> BMP-2 stimulates calcification (red spots, left column), whereas cells cultured in the presence of PODS<sup>®</sup> Empty crystals or cells alone did not show signs of osteogenic differentiation. Only the cells in close vicinity of the PODS<sup>®</sup> BMP-2 field stained red, demonstrating build-up of BMP-2 around the PODS<sup>®</sup> crystals.

**B**

PODS <sup>®</sup> BMP-2	PODS <sup>®</sup> Empty	Cells only
----------------------------	----------------------------	---------------



**(B) Alcian blue staining of ATDC5 cells**

ATDC5 cells stained with Alcian blue to indicate chondrogenic differentiation (glycosaminoglycan synthesis). Analogous to the Alizarin red staining, only PODS<sup>®</sup> BMP-2 crystals were able to stimulate ECM production in ATDC5 cells (blue spots, left column), whereas cells that were cultured in the presence of PODS<sup>®</sup> Empty crystals or cells alone did not show signs of differentiation towards cartilage formation.

## Conclusions

- PODS<sup>®</sup> crystals adhere efficiently to plastic surfaces, ideal for localized placement on tissue culture dishes.
- PODS<sup>®</sup> achieve growth factor gradients easily when positioned in high densities in confined areas.
- Even for long culture periods, a single application of PODS<sup>®</sup> crystals is effective, significantly reducing both hands-on time and cost of materials.
- PODS<sup>®</sup> BMP-2 exhibits bioactivity.

For more information and a full list of our current PODS<sup>®</sup> growth factors, please visit our website [www.cellgs.com](http://www.cellgs.com).



---

Cell Guidance Systems' reagents and services enable control, manipulation and monitoring of the cell, both *in vitro* and *in vivo*

**Growth Factors**

- Recombinant
- Sustained Release

**Exosomes**

- Purification
- Detection
- Tracking
- NTA Service

**Small Molecules**

**Cell Counting Reagent**

**Matrix Proteins**

**Cell Culture Media**

- Pluripotent Stem Cells
- Photostable
- *In Vitro* Blastocyst Culture
- ETS-embryo Culture
- Custom Manufacturing Service

**Gene Knock-Up System**

**Cytogenetics Analysis**



General [info@cellgs.com](mailto:info@cellgs.com)  
Technical Enquiries [tech@cellgs.com](mailto:tech@cellgs.com)  
Quotes [quotes@cellgs.com](mailto:quotes@cellgs.com)  
Orders [order@cellgs.com](mailto:order@cellgs.com)

[www.cellgs.com](http://www.cellgs.com)

**EUROPE**  
**Cell Guidance Systems Ltd**  
Maia Building  
Babraham Bioscience Campus  
Cambridge  
CB22 3AT  
United Kingdom  
T +44 (0) 1223 967316  
F +44 (0) 1223 750186

**USA**  
**Cell Guidance Systems LLC**  
Helix Center  
1100 Corporate Square Drive  
St. Louis  
MO 63132  
USA  
T 760 450 4304  
F 314 485 5424