



# Competent Cells 101:

Maximizing Your Success in Cloning, Protein  
Expression, and Library Generation

# Agenda

- Selecting your competent cell format:
  - electrocompetent vs. chemically competent
- Maximizing transformation efficiency to create large libraries
- Choosing the right strain:
  - Strain genotypes and properties
  - Specialized strains for specific applications
- Engineered strains for low-endotoxin protein and plasmid preparation



# Competent Cell Choices

Multiple differences between formats

## Chemically Competent Cells:

- Easy to work with, less expensive
- Do not require specialized equipment
- Accept larger DNA volumes
- Lower transformation efficiency than Electrocompetent Cells
- Used for routine cloning and plasmid propagation

## Electrocompetent Cells:

- Shorter protocol, amenable to automation
- Highest transformation efficiency
- Require an electroporator instrument
- Require low DNA volumes
- Used for library construction and propagation




# How Does Transformation Work?

Protocols are optimized for maximum cell viability and efficiency

## Chemically Competent (Heat Shock) Transformation:

- Cells and DNA are incubated together on ice
- Salt masks the negative charge on DNA
- Heat shock (42°C) enhances DNA uptake
- Media is added and cells are “recovered” for 1 hour at 37°C



## Electrocompetent Transformation:

- Cells and DNA are mixed immediately before electroporation
- Electrical current is applied, creating pores in the membrane
- Cold temperature is maintained to reduce cell damage
- Media is added and cells are “recovered” for 1 hour at 37°C



# Transformation Workflow Comparison

Protocols have parallel steps, different time and equipment requirements

## Chemically Competent Cells

1. Remove cells from  $-80^{\circ}\text{C}$  and thaw on ice
2. Add DNA to cells
3. Incubate on ice 30 min
4. Heat shock  $42^{\circ}\text{C}$  45 sec.
5. Incubate on ice 2 min
6. Add Recovery Medium
7. Incubate 1 hour at  $37^{\circ}\text{C}$  with shaking



~ 1 hour 45 minutes

## Electrocompetent Cells

1. Remove cells from  $-80^{\circ}\text{C}$  and thaw on ice
2. Add DNA to cells
3. Transfer cells + DNA to pre-chilled electroporation cuvette
4. Perform electroporation using recommended parameters
5. Immediately add Recovery Medium
6. Incubate 1 hour at  $37^{\circ}\text{C}$  with shaking




~ 1 hour 10 minutes



# What is Transformation Efficiency?

Highest transformation efficiency is achieved with very low amounts of DNA

- Measures the rate at which cells take up plasmid and form colonies
- Reported as Colony Forming Units per  $\mu\text{g}$  DNA (CFU/ $\mu\text{g}$ )
- Determined using intact, supercoiled plasmid DNA (10 pg of pUC19 plasmid is standard)



In practice, higher [DNA] is used for library construction






# Why is Efficiency Important?

*Higher efficiency = Higher Diversity = Better Library Screen = More Confidence in Results*

DNA libraries are **collections of DNA clones**

Biotech/pharma customers **require high efficiencies** for library construction


Customers need the highest efficiency cells possible to produce the **highest diversity libraries**



By Aluquette (Own work) [CC BY-SA 3.0 (<http://creativecommons.org/licenses/by-sa/3.0>)], via Wikimedia Commons



# Tips for Library Transformations



- Direct electroporation of ligation reactions
  - *heat inactivate ligase as recommended by supplier*
- Clean up ligation reaction to remove ligase and buffer
- Recover DNA in water at 200 ng/ $\mu$ L or higher
- Use up to 200-300 ng ligated DNA in 1  $\mu$ L H<sub>2</sub>O per standard 25  $\mu$ L electroporation reaction
- Scale-up: maintain 200-300 ng DNA in 1  $\mu$ L H<sub>2</sub>O per 25  $\mu$ L cells
  - *test and optimize at small scale prior to scale-up*
- For cell volumes 50  $\mu$ L to 300  $\mu$ L, use 2 mm cuvettes
  - *optimization of electroporation parameters may be needed*





# Plating Medium Matters

Salt concentration is *critical* for optimal colony growth.


*All LB Medium is not created equal.*

## LB-Lennox (per liter)

5 g NaCl, Low Salt

5 g Yeast Extract

10 g Tryptone




## LB-Miller (per liter)

10 g NaCl, High Salt

5 g Yeast Extract

10 g Tryptone



*E. Cloni*<sup>®</sup> 10G cells overnight growth comparison (~16 hours).

Use **Low-Salt LB-Lennox** for most Lucigen cloning and protein expression strains

(\*Only exception: ClearColi<sup>®</sup> competent cells should be plated using *High-Salt* LB-Miller).



# Competent Cell Packaging Formats

Choose the best solution for your workflow

## Routine cloning and library construction:


- 1 reaction per tube (SOLO)
- 2 reactions per tube (DUO)
- 6 reactions per tube (6-Packs)

## High-throughput cloning:

- 96-well microarray plates

## Custom formats for any workflow:

- Customer-specified vessel formats, aliquot sizes and strains





**What competent cell strain(s)  
do you use in your lab?**



# Lucigen Competent Cell Categories

Cells for all applications

- General cloning and library construction
- Phage Display library construction
- Large/Difficult fragment cloning
- Protein Expression
- Low Endotoxin



# General cloning and library construction

Strain	Supports cloning Methylated DNA	Contains F' Plasmid	Suitable for Blue/White Screening	IPTG needed for Blue/White Screening
<i>E. coli</i> ® 10G	Yes	No	Yes	No
<i>E. coli</i> 10GF'	Yes	Yes	Yes	Yes
<i>E. coli</i> 5-alpha	No	No	Yes	No



**Mammalian DNA**  
**Plant DNA**



# Genotypes of General Cloning Strains

***E. cloni*® 10G:** F- ***mcrA*Δ(*mrr-hsdRMS-mcrBC*) *endA1 recA1* Φ80 Δ(*lacZ*)M15**  
 Δ*lacX74 araD139* Δ(*ara,leu*)7697 *galU galK rpsL nupG* λ- ***tonA*** (StrR)

***E. cloni* 10GF':** [**F'** *proA+B+* ***lacI<sup>q</sup>Z*ΔM15::Tn10 (TetR)**] / ***mcrA* Δ(*mrr-hsdRMS-***  
***mcrBC*) *endA1 recA1* Φ80 Δ(*lacZ*)M15** Δ*lacX74 araD139*  
 Δ(*ara,leu*)7697 *galU galK rpsL nupG* λ- ***tonA*** (StrR)

***E. cloni* 5-alpha:** ***fhuA2*Δ(*argF-lacZ*)U169 *phoA glnV44* Φ80 Δ(*lacZ*)M15 *gyrA96***  
***recA1 relA1 endA1 thi-1 hsdR17***


Genotype	Function
<i>endA1</i>	Eliminates endonuclease that commonly contaminates DNA preps and reduces yield
<i>recA1</i>	Reduced recombination
<i>tonA, fhuA2</i> Δ	Resistance to T1 bacteriophage
F'	Allows M13 infection for ssDNA production; Carries <i>lacI<sup>q</sup></i> , overexpresses lac repressor
<i>mcrA</i> Δ( <i>mrr-hsdRMS-mcrBC</i> )	Eliminate restriction of methylated foreign DNA
Φ80 Δ( <i>lacZ</i> )M15	β-galactosidase ω fragment – blue-white screening





# Phage Display Strains

- SS320
- TG1
- ER2738
- MC1061 F-



# Phage Display Strains


Cell Strain	Uses	Amber Suppressor	Transformation Efficiency
TG1	Phage Display Library Screening Protein Expression	Yes	$\geq 4 \times 10^{10}$ cfu/ $\mu$ g
SS320 (MC1061F')	Phage Display Library Screening	No	$\geq 4 \times 10^{10}$ cfu/ $\mu$ g
MC1061 F- (same as SS320 but lacks F' episome)	General Cloning Phage Display	No	$\geq 4 \times 10^{10}$ cfu/ $\mu$ g
ER2738	Phage Display Library Screening	Yes	$\geq 2 \times 10^{10}$ cfu/ $\mu$ g



# Cloning Large/Difficult DNA

## What makes DNA difficult to clone?

- Toxic coding sequences
- Promoters
- A-T Rich DNA
- Large fragments (>10 kb)
- Repetitive sequences




## What are the solutions?

- Specialized competent cell strains (e.g. Endura)
- Vectors and Vector/Host systems (e.g. BigEasy)

# Endura™ Competent Cells

Clone unstable sequences

- Ideal for construction of lentiviral vector libraries
- Reduce unwanted recombination
- Stabilize DNA repeats
- Obtain high DNA yields
- Direct replacements for Stbl3™ and SURE strains
- Available in electrocompetent or chemically competent formats



# Endura™ Competent Cells

Recommended by the Broad Institute for  
lentiviral GeCKO CRISPR Libraries



LENTIVIRAL CRISPR TOOLBOX




## Genome-scale CRISPR Knock-Out (GeCKO) v2.0 pooled libraries

CRISPR (Clustered Regularly Interspaced Short Palindromic Repeats) is a microbial nuclease system containing a combination of CRISPR-associated (Cas) genes as well as non-coding RNA elements capable of programming the specificity of the CRISPR-mediated nucleic acid cleavage. Using lentivirus, we delivered the type II CRISPR nuclease system to facilitate genome editing in mammalian cells in a pooled library targeting early consecutive exons ([Shalem\\*, Sanjana\\*, et al., Science 2014](#)). Here we describe how to amplify GeCKO v2.0 DNA plasmids to have sufficient quantity to produce lentivirus, while maintaining full library representation.

**Library description:** The GeCKO v2 libraries consist of specific sgRNA sequences for gene knock-out in either the human or mouse genome. Each species-specific library is delivered as two half-libraries (A and B). When used together, the A and B libraries contain 6 sgRNAs per gene (3 sgRNAs in each library). We recommend screening the entire library (A and B) when possible but if adequate representation cannot be obtained with the entire library, screening can be performed with one half-library. Since more cells are needed in the screen as number of constructs in the library increases, this design has the flexibility of screening with just the A library (3 sgRNAs per gene) at a smaller scale or screening the full library (6 sgRNAs per gene). Both A and B libraries contain 1000 control sgRNAs designed not to target in the genome. The A library also targets miRNAs (4 sgRNAs per miRNA). Each library is available in a 1 vector (lentiCRISPRv2) or 2 vector (lentiCas9-Blast and lentiGuide-Puro) format. The 2 vector format with the library in lentiGuide-Puro has the advantage of higher titer for the library virus but requires cells to already contain Cas9 (usually genomically integrated using lentiCas9-Blast). See p. 3 for full library specs.


### One vector lentiviral GeCKO system



### Two vector lentiviral GeCKO system



Visit: <http://bit.ly/CRISPR-Endura>



# Bacterial Vector/Host Systems

## Systems for Any Insert Size



### Options available:

- CloneSmart® kits for <15 kb
- BigEasy® v2.0 kits for exceptionally difficult targets
- CopyRight® kits for Fosmid or BAC cloning

### *Convenient, fast, easy, and effective*

- Pre-cut, dephosphorylated – ready to clone
- <1% Background
- Minimize cloning problems

See recent webinar on “Cloning the Unclonable”

<http://www.lucigen.com/webinars.html>






# Lucigen Competent Cells for Protein Expression

## T7 Expression hosts

- *E. coli*<sup>®</sup> EXPRESS BL21(DE3): Electrocomp and Chem Comp
- HI-Control™ BL21(DE3) Chemically Competent Cells
- OverExpress™ C41(DE3) & C43(DE3): Electrocomp and Chem Comp
- ClearColi BL21(DE3): Electrocompetent Cells

## Non-T7 hosts for expression from bacterial promoters

- *E. coli* 10G
- HI-Control™ 10G
- Biotin XCell™ F'



# Cloning and Expression of Toxic Targets

Background expression of toxic targets is problematic in cloning and expression

## Problem:

Background expression of T7 RNA Polymerase

= *Uninduced target gene expression*


*Decreased colony count and cell viability*

## Solutions:

HI-Control BL21(DE3)

OverExpress C41(DE3) and C43(DE3)


*Improved control of T7 RNA polymerase*



# OverExpress<sup>®</sup> C41(DE3) & C43(DE3)

Better Viability and Expression of Toxic & Membrane Proteins

- BL21(DE3) derivatives
- Useful for expression from T7 promoters
- Cited in over 350 publications



Dumon-Seignovert, et al. , (2004). *Protein Expression and Purification* **37**, 203-206.



# Competent Cells for Low Endotoxin

Protein and Plasmid Production


Do you currently need  
endotoxin-free DNA (plasmid)  
and/or protein preps  
in your research?



# What is Endotoxin Contamination?

Endotoxin is present in E. coli protein and DNA preps


- **Endotoxin, or LPS (lipopolysaccharide)**, is part of the bacterial outer membrane
- LPS is problematic:
  - Elicits an immune response in human cells
  - Affects cell viability and growth
- LPS **must** be removed from therapeutic proteins, cell-based assay ligands and cell culture additives
- **Endotoxin removal is difficult, costly and time-consuming**



# ClearColi<sup>®</sup> Does Not Produce LPS

Save time and money by eliminating downstream LPS removal

- ClearColi is genetically engineered to produce Lipid IV<sub>A</sub> instead of LPS
- Contains 7 deletion mutations for modifying LPS to Lipid IV<sub>A</sub> and 1 compensating mutation
- Structure of Lipid IV<sub>A</sub> *does not* trigger an immune response in human cells






# HEK-Blue™ hTLR4/MD-2 Stimulation Assay

Lipid IV<sub>A</sub> does not trigger the endotoxic response in human cells

## LPS vs. Lipid IV<sub>A</sub>





No Stimulation  
from Lipid IV<sub>A</sub>



# ClearColi<sup>®</sup> BL21(DE3) Strain for Protein Production

- Achieve equivalent protein expression yield & profile
- Reduce endotoxin by > 97%
- Lipid IV<sub>A</sub> can be removed during standard purification (ie, Ni<sup>2+</sup> column)



# ClearColi<sup>®</sup> K12 Strain for Plasmid Production

- Produce low-endotoxin plasmid DNA
- Achieve equal or better plasmid DNA yield
- Reduce endotoxin by >99% using standard plasmid prep kits

	Plasmid Yield per 1 mL
ClearColi K-12	4.83 µg
<i>E. coli</i> <sup>®</sup> 10G	3.75 µg

Endotoxin Levels Detected  
ClearColi<sup>®</sup> K-12 vs. DH10B *E. coli*  
EU/mg DNA



# Important ClearColi Facts

- ClearColi cells require higher NaCl than other Lucigen strains
  - use LB Miller medium (10 g/L NaCl)
- ClearColi cells grow more slowly than standard strains
- LAL assay has some sensitivity to lipid IV<sub>A</sub>, but human TLR4/MD2 does not
- Different mammalian systems have different sensitivities to lipid IV<sub>A</sub>

<http://www.lucigen.com/faq-clearcoli.html>



# Custom Competent Cells

**Your strain of choice, made to the highest competency**

- ✓ Made-to-order competent cells
- ✓ Lucigen's strains and customer-supplied strains
- ✓ Chemically competent or electrocompetent cells
- ✓ Custom dispense volumes and formats
- ✓ Customer-specified tubes and plate types
- ✓ Up to  $4 \times 10^{10}$  cfu/ $\mu$ g
- ✓ Lucigen is registered to the ISO13485 quality management standard for medical devices through BSI, Inc.





# Sample Lucigen Competent Cells

**Free samples for the following strains:**

*E. cloni*<sup>®</sup> 10G ElectroComp Cells

*E. cloni* 10G ChemComp Cells

Endura<sup>™</sup> ElectroComp

Endura ChemComp

Phage Display ElectroCombo Pack

Includes TG1, SS320 (MC1061 F') and ER2738

<http://lucigen.com/samplecells.html>



eLucidator<sup>™</sup>  
**Competent Cells**  
INTERACTIVE SELECTION GUIDE



Find the right competent cells for your application:

<http://lucigen.com/compcellguide>





Questions? [www.lucigen.com](http://www.lucigen.com)

Thank You!

Lucigen Tech Support  
[techsupport@lucigen.com](mailto:techsupport@lucigen.com)  
(888) 575-9695  
(608) 831-9011  
8 am – 5 pm Central Time

