

PROTOCOL FOR PJAZZ® CLONES



METHOD FOR HIGH THROUGHPUT SEQUENCING OF PJAZZ CLONES

Protocol courtesy of Laboratory for Genomics and Bioinformatics, University of Oklahoma Health Sciences Center, OK
<http://www.microgen.ouhsc.edu>

1 Prep DNA using Millipore Montage Plasmid Miniprep 96 well kit.

2 Set up the sequencing reactions using Phenix 384-well FrameStar PCR plates.

- **On ice, make a 2X Master Mix for 100 rxns:**
 - 16.50 μ l Big Dye V3.1
 - 0.78 μ l Primer @ 320 pmol/ μ l (2.5 pmol/rxn)
 - 118.00 μ l 5X ABI sequencing buffer
 - 164.72 μ l DIUF water
 - 300.00 μ l Total
 - **Use 3 μ l of 2X Master Mix per rxn**
 - **Add 3 μ l of template (60ng-180ng) per rxn**
-

3 Cycle.

- **95 °C for 4 min**
 - **then 25 cycles of:**
 - 95 °C for 15 sec
 - 55 °C for 5 sec
 - 60 °C for 2 min
-

4 Hold at 4 °C until ready for use.

5 Clean up with ethanol precipitation or Sephadex (G-50).
