


Separate



**SurePage Pre-cast Gels**

Fast run times of 20+ Min  
High performance  
Free running buffer



**Protein Standards**

Comprehensive selection & Broad  
Molecular Weight Range

Transfer



**eBlot™ L1 Protein Transfer System**

Achieve complete transfer  
in 9-17 Min



**eStain™ L1 Protein Staining System**

Ponceau staining  
in 10 Min

Wash  
&  
Incubate



**eZwest™ Lite Automated Western Device**

Automated Western Blotting device completing the entire  
process of blocking, antibody incubation and washing steps.

Image  
&  
Quantify



**THE™ Elite Antibodies**

Precise and effective antibodies that meet the most  
demanding research requirements.

# ONE-STOP SOLUTION FOR WESTERN BLOTTING

## Ordering Information

### SurePAGE Pre-Cast Gels

Cat No.	Product Name	Quantity
M00661	SurePAGE™, Bis-Tris, 10×8, 8%, 10 wells	10 PK
M00662	SurePAGE™, Bis-Tris, 10×8, 8%, 12 wells	10 PK
M00663	SurePAGE™, Bis-Tris, 10×8, 8%, 15 wells	10 PK
M00664	SurePAGE™, Bis-Tris, 10×8, 10%, 10 wells	10 PK
M00665	SurePAGE™, Bis-Tris, 10×8, 10%, 12 wells	10 PK
M00666	SurePAGE™, Bis-Tris, 10×8, 10%, 15 wells	10 PK
M00667	SurePAGE™, Bis-Tris, 10×8, 12%, 10 wells	10 PK
M00668	SurePAGE™, Bis-Tris, 10×8, 12%, 12 wells	10 PK
M00669	SurePAGE™, Bis-Tris, 10×8, 12%, 15 wells	10 PK
M00652	SurePAGE™, Bis-Tris, 10×8, 4-12%, 10 wells	10 PK
M00653	SurePAGE™, Bis-Tris, 10×8, 4-12%, 12 wells	10 PK
M00654	SurePAGE™, Bis-Tris, 10×8, 4-12%, 15 wells	10 PK
M00655	SurePAGE™, Bis-Tris, 10×8, 4-20%, 10 wells	10 PK
M00656	SurePAGE™, Bis-Tris, 10×8, 4-20%, 12 wells	10 PK
M00657	SurePAGE™, Bis-Tris, 10×8, 4-20%, 15 wells	10 PK
M00658	SurePAGE™, Bis-Tris, 10×8, 8-16%, 10 wells	10 PK
M00659	SurePAGE™, Bis-Tris, 10×8, 8-16%, 12 wells	10 PK
M00660	SurePAGE™, Bis-Tris, 10×8, 8-16%, 15 wells	10 PK

### Protein Standards

Cat No.	Product Name	Quantity
M00624-250	Broad Multi Color Pre-Stained Protein Standard	250 µL
M00624-1250	Broad Multi Color Pre-Stained Protein Standard	1250 µL
M00516	PAGE-MASTER Protein Standard	500 µL
MM1397-500	PAGE-MASTER Protein Standard Plus	500 µL
M00521	WB-Master Protein Standard	250 µL

### eStain L1 Protein Staining System

Cat No.	Product Name	Quantity
L00657	eStain L1 Protein Staining Device	1 Unit
L00753	eStain L1C Protein Staining Kit	1 Kit
L00659-1	eStain L1 Protein Staining Kit	1 Kit

### eBlot L1 Fast Transfer System

Cat No.	Product Name	Quantity
L00686	eBlot L1 Fast Transfer Device	1 Unit
L00724	NC eBlot L1 Transfer Kit (Excludes NC Membranes)	1 Kit
L00725	Complete eBlot L1 Transfer Kit (NC)	1 Kit
L00726	PVDF eBlot Transfer Kit (Excludes PVDF Membranes)	1 Kit
L00727	Complete eBlot L1 Transfer Kit (PVDF)	1 Kit

### eZwest Automated Western Device

Cat No.	Product Name	Quantity
L00816	eZwest Automated Western Device	1 Unit

For more information, visit us at [www.genscript.com](http://www.genscript.com)