

Steroid Hormone Assay Kits



Researchers have trusted Cayman's collection of steroid hormone ELISAs for decades to study effects in metabolism, inflammation, immune function, osmoregulation, sexual development, and stress responses. Designed with competitive assay principals, each kit offers all reagents needed for specific quantification of steroids in cell culture medium, plasma, and serum.

Assay Features:

- Highly sensitive detection in the low pg/ml range
- Validated in several biological matrices
- Broad species reactivity

Androgens

Item No.	Product Name
582701	Testosterone ELISA Kit
582751	11-keto Testosterone ELISA Kit
501740	Testosterone Glucuronide ELISA Kit

Estrogens

Item No.	Product Name
501890	Estradiol ELISA Kit
582281	Estriol ELISA Kit
501290	Estrone-3-Glucuronide ELISA Kit
582261	2-Methoxyestradiol ELISA Kit

Progestins

Item No.	Product Name
501300	Pregnanediol-3-Glucuronide (PDG) ELISA Kit
582601	Progesterone ELISA Kit

Glucocorticoids

Item No.	Product Name
501320	Corticosterone ELISA Kit
500360	Cortisol ELISA Kit
500370	Cortisol Express ELISA Kit
501420	11-Oxoetiocholanolone ELISA Kit

Mineralocorticoids

Item No.	Product Name
501090	Aldosterone ELISA Kit

Vitamin D

Item No.	Product Name
501050	Vitamin D ELISA Kit

Newly Updated Estradiol ELISA Kit - Item No. 501890

Our scientists developed a new antiserum with less cross reactivity compared to the previous kit and changed the sample preparation/extraction protocol to a simpler and less hazardous format that provides improved specificity.

View a side-by-side comparison of the previous and new kit at www.caymanchem.com/estradiol

Find answers to frequently asked questions about our steroid hormone kits at www.caymanchem.com/steroidFAQ